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Rev. n°: 11/2016

PREVENTS THE OCULAR DETERIORATION ASSOCIATED TO AGING AND HARSH ENVIRONMENTAL CONDITIONS,
CONTRIBUTING TO THE MAINTENANCE OF VISUAL ACUTY

- Clear and sharp vision
- Vasoprotector

- Antioxidant
- Moisturiser

Preservative free
Sterile solution

Preservative-free solution which, thanks to its high technology filter, remains sterile throughout its continued use. As it contains no preservatives, it is gentle on the surface of the eye and can be used frequently.

visglyc® neo is a sterile ophthalmic solution, with a physiological pH, based on natural substances: N-acetylcarnosine, *Vaccinium myrtillus* extract and chondroitin sulphate of marine origin.

The combined action of these natural substances:

- Prevents ocular damage derived from aging and harsh environmental conditions contributing to the maintenance of the visual acuity, with antioxidants such as N-acetylcarnosine and *Vaccinium myrtillus* extract that aid in scavenging free radicals that damage the ocular crystalline lens.
- Contributes to regulate the microcirculation and ocular nutrition, due to the hydration capacity of chondroitin sulphate and the vasoprotector property of *Vaccinium myrtillus*.

INSTRUCTIONS FOR USE

Open the bottle and squeeze gently to apply 1 drop in each eye. Replace cap after each use. Since it contains no preservatives, it can be used frequently. Compatible with contact lenses. Can be use up to 3 months after opening. Container for personal use.

COMPOSITION

N-acetylcarnosine (1%), *Vaccinium myrtillus* extract, chondroitin sulphate (0.15%), sodium phosphate dibasic dodecahydrate, sodium phosphate monobasic monohydrate, sodium chloride, water for injections q.s. 100ml.

PRECAUTIONS AND WARNINGS

- Exclusively for ophthalmic use.
- Wash your hands before application.
- Do not use the product if the container, before opening, is damaged or not well closed.
- Do not use in case of sensitivity or allergy to any of the components.
- During application, it is recommended not to touch the eye or any other surface with the container applicator.
- Do not require any special storage conditions.
- Do not use after the expiry date indicated on the container.
- Keep out of reach of children.
- Do not dispose of the container in the environment after use.
- It is advised to consult a doctor about the compatibility with other medicines.
- If discomfort is noted, suspend the treatment and consult a specialist.

PRESENTATION

Multidose container with 10 ml of sterile solution.

SHELF LIFE

Unopened, do not use after the expiry date stated on the box and the bottle label. The contents can be used up to 3 months after opening.

Rev nº: 11/2016

MEDICAL DEVICE

STERILE | A

MANUFACTURER:
OPKO Health Spain S.L.U.
Pl. Europa 13-15, Local 2,
08908 L'Hospitalet de Llobregat
(Barcelona), Spain